

Understanding Protective stops

*Applicable to robots running SW 5.6 or
newer*

Oluf Skov Rosenlund, Tech Lead, Robotics,
Product Creation, Universal Robots

Document version 1.0 24/10 2023





Understanding Protective Stops

This presentation will hopefully answer some of these questions:

What exactly is a protective stop?

What are typical application problems that result in protective stops?

When the system fails, what does the errors tell you?

How to analyze and solve protective stop problems?

Agenda

1. Faults, Violations and Protective stops
2. A bit of tech
3. Protective stops
 - Protective stops due to robot (correctly or incorrectly) detects a collision
 - Protective stops that stops the robot before target trajectory exceeds a safety limit
 - Other protective stops
4. Questions?





Faults, Violations and Protective stops

Protective stops

What it is:

A protective stop is issued when the robot can no longer perform the intended motion.

Detection:

Controller detects that the target trajectory cannot be followed.

Effect:

Robot decelerates, program can be resumed.

Note:

Different types of protective stop have different root causes and mitigations.

Robot can be freedriven or just pushed if it has collided with something.



Safety system Violations

Detection:

Safety system observes a safety limit is exceeded.

Effect:


Robot decelerates and powers off motors, program can be resumed.

Note:

This should in theory never happen, Controller ensures to always react before safety system activates.

(the error on the picture usually happens if you pull the robot hard while it is e.g. Safeguard stopped).

The robot should just be restarted and the application can resume.



Safety system Faults

Detection:

Safety system detects a fault in hardware or software.

Effect:

Robot decelerates and powers off motors.


Note:

Faults can be caused by incorrectly wired safety IO (causes safety system disagreements)

Try installing latest software and see if problem is resolved (in case of a software issue)

If problem occurs multiple times, component may need repair

See log in Log Viewer for details




Safety Stop approach

Protective stops, Safeguard stops, E-stops and violations all pause the program so it can be resumed from where it was paused.

Robot always stops within the *stopping time* and *stopping distance* safety limits

Limit	Normal	Reduced
Power	1000 W	300 W
Momentum	100.0 kg m/s	25.0 kg m/s
Stopping Time	1000 ms	400 ms
Stopping Distance	2000 mm	500 mm
Tool Speed	5000 mm/s	1500 mm/s
Tool Force	250.0 N	150.0 N
Elbow Speed	5000 mm/s	1500 mm/s
Elbow Force	250.0 N	150.0 N



Flow of E-stops, Violations and Faults (Cat. 0 and 1 stops)*

1. Decelerate movement on trajectory down to zero velocity
2. Engage brakes when standing still
3. Remove power from motors so robot "parks" on the brake (robot drops a few mm during that process)

Flow of Protective stops, Safeguard stops and 3 position enabling device (Cat. 2 stops)*

1. Decelerate movement on trajectory down to zero velocity
2. If a collision is detected (C153, C157, C158, C159), robots performs a small reverse movement

*except very rare fault scenarios where encoders or motor is broken and a Cat.0 stop

A bit of tech

Tech: Torque windows

- Force and power limits affect allowed torque/current in motors
- Window is affected by robot pose
- Torque window is centered around target current




Limit	Normal	Reduced
Power	1000 W	300 W
Momentum	100.0 kg m/s	25.0 kg m/s
Stopping Time	1000 ms	400 ms
Stopping Distance	2000 mm	500 mm
Tool Speed	5000 mm/s	1500 mm/s
Tool Force	250.0 N	150.0 N
Elbow Speed	5000 mm/s	1500 mm/s
Elbow Force	250.0 N	150.0 N

Force Limits Example

Example base torque limits:
Tool and elbow force limit = 150 N

Base-tool dist: 0.9 m
 $0.9\text{m} \times 150\text{N} = 135\text{ Nm}$
 Base-elbow dist: 0.5 m
 $0.5\text{m} \times 150\text{N} = 75\text{ Nm}$

Base torque window: 75Nm





Tech: Time scaling

MoveJ and MoveL in SW5.0 and newer:

Same velocity and acceleration for two moves.

Movement acceleration is automatically scaled due to payload change



Which ways do we have to know if safety system is taking effect on any moment of the cycle? Is there a way to identify which of the safety parameters is actually affecting?

It is quite hard to see what limit is scaling the robot.

- During accelerations/decelerations it is fixed limits the prevent exceeding mechanics and gear limits, or PSU limits (e.g. for the OEM CB)
- During cruise phases, it can be joint speed, TCP speed, elbow speed, momentum.
- During vertical movements typically, it is stopping time or stopping distance (a stretched robot with full payload will move slowly downwards as it almost does not have enough torque to brake if needed)



Protective stops

Types of protective stops

A protective stop is issued when the robot can no longer perform the intended motion.

Different types of protective stop have different root causes and mitigations.

Protective stops due to robot (correctly or incorrectly) detects a collision:

- C153: Position deviates from path
- C159: Position deviates from path (zero payload)
- C157: Collision detected by joint
- C158: Collision detected by joint (zero payload)

Protective stops that stops the robot before target trajectory exceeds a safety limit:

- C150: Position close to joint limits
- C151: Tool orientation close to limits
- C152: Position close to safety plane limits

Other protective stops:

- C154: Position in singularity
- C155: Robot cannot maintain its position, check if payload is correct
- C156: Wrong payload or mounting detected, or something is pushing the robot when entering Freedrive mode

Protective stops due to robot
(correctly or incorrectly)
detects a collision

Error codes:

C153: Position deviates from path

C159: Position deviates from path
(zero payload)

C157: Collision detected by joint

C158: Collision detected by joint
(zero payload)

Check:

Deviation between target and actual
position and velocity

Note:

Root causes:

- Wrong payload, wrong mounting
- External cabling dragging the robot, or collisions
- Flexible mounting or vibrating tools
- Robot on external axis
- Very high accelerations combined with blends
- External control with `servoj()` movements with very high accelerations and decelerations, or jitter
- Startup in cold environments

Let's dig into what you can do about those root causes

Protective stops due to wrong payload or mounting

Detection:

Test if the payload is correct by setting the robot in freedrive. If the robot feels gravity compensated in multiple positions and orientation, these settings are good.

Note:

TCP offset does **not** affect the dynamic model of the robot, so it does not influence protective stops.

Guidelines:

- Set payload, Center of Gravity in installation settings
- If picking and placing objects heavier than 1kg, payload should be adjusted during program using Payload program node in PolyScope or `set_target_payload(m, CoG)` script command
- For picking heavy objects with suction cup grippers you may use **payload transition time** to smoothly adjust the payload
`set_target_payload(m, cog, inertia=[0, 0, 0, 0, 0, 0], transition_time=0)`


See also <https://www.universal-robots.com/articles/ur/application-installation/how-to-design-for-higher-payloads/>

Protective stops due to collisions

Detection:


Actual current suddenly deviates from target.

Position error builds up fast



Guidelines:

- Visually inspect what caused the collision
- Use "move until contact" PolyScope node or Force control when e.g. stacking.




- Use urscript command `position_deviation_warning()` to check if parts of a program is close to causing a protective stop
- Higher force limits for tool and elbow increases allowed motor torque to move on the trajectory

Protective stops due to vibrations


Detection:

Robot vibrates, oscillations in realltime data



Guidelines:

- Ensure stiff mounting
- Set payload inertia in PolyScope




See also <https://www.universal-robots.com/articles/ur/application-installation/how-to-design-for-higher-payloads/>

Protective stops due to accelerations of external axes

Detection:

Varying current consumption when external axis accelerates



Guidelines:

New in SW 5.14:

Use `set_base_acceleration(a)` script command to compensate for acceleration

Use `high_holding_torque_disable()` to make the robot detect collisions when the axis is moving and the robot is standing still

In older software, you can use `set_gravity()` script command

See also new commands in <https://www.universal-robots.com/articles/ur/programming/set-base-acceleration-and-high-holding-torque-disableenable/>

Protective stops due to very high accelerations combined with blends

Notes:

MoveJ and MoveL without blends scales torques to be within hardware limits


Detection:

Exceeding max torque during blends
Look for warning C173 joint torque overload

Guidelines:

- Adjust accelerations to lower required torque during blends
- Scaled MoveJ and MoveL gives **optimal** cycle times
- In SW 5.6 and newer, robot will do a controlled deceleration within torque limits if leaving a “continuous if” statement or having a blend in the last waypoint of a move command. It is still a good idea to add a `stopj()` or a `stopl()` if you want a lower acceleration.

Look out for future software releases with improved handling of blends



External control with `servoj()` movements with very high accelerations and decelerations, or jitter

Notes:

- * Accelerations are scaled to be within torque limits, and stay on trajectory
- * Decelerations are the challenge

Detection:

Robot overshooting waypoints, or too high torque spikes during decelerations, or protective stops during decelerations



Guidelines:

- Consider if you can switch to MoveJ or MoveL, then robot controller scales accelerations optimally for you.
- If possible, run real-time part of trajectory generation on robot controller, either as a daemon or in urscript to avoid jitter. Then you can feed positions with 500Hz
- Adjust decelerations in your trajectory generator software
- Drivers exist here:
 - https://github.com/UniversalRobots/Universal_Robots_ROS_Driver
 - Or https://gitlab.com/sdurobotics/ur_rtde


These drivers do not resolve the challenge of choosing optimal accelerations yet.

Drivers can be used as inspiration; extend with the functions you need!

Startup in cold environments

Detection:

Slow buildup of position and velocity error



Guidelines:

- With warmup program, robot can work down to minus 15°C ambient
- Make program with slow warmup movement
- Use SW 5.13.1 or newer
- Tested in climate chamber at UR



Protective stops that stops the robot before target trajectory exceeds a safety limit

Error codes:

C150: Position close to joint limits

C151: Tool orientation close to limits

C152: Position close to safety plane limits

Root cause:


Target trajectory exceeds a safety limit.

Robot software prevents exceeding the limit by issuing a protective stop.

Can happen if safety planes are enabled/disabled in reduced mode.

Guidelines:

Check that waypoints are placed on the right side of safety planes



Other protective stops: C155 and C156: Robot cannot maintain its position

C155: Robot cannot maintain its position, check if payload is correct

Detection:

Happens if a very wrong payload or mounting is entered. Robot starts moving upwards or downwards due to wrong dynamics. Movement is detected and stop is issued

Guideline:

Make sure payload and mounting is correct

Is something external pushing the robot?

C156: Wrong payload or mounting detected, or something is pushing the robot when entering Freedrive mode

Detection:

The robot software checks if a sudden acceleration happens right after entering freedrive. If that happens, it protective stops as this could be caused by wrong payload setting

Guideline

The check is a bit sensitive...

Try to avoid pulling the robot before the freedrive button is pressed

Don't use freedrive when mounting or unmounting tools, unless you know what you do...

C154: Position in singularity

C154: Position in singularity

Detection:

Happens if the target trajectory + offsets from conveyor tracking, force mode, or path offset

Guidelines:


Consider if base should be tilted to avoid base singularity

Consider mounting tool at an angle to avoid wrist alignment singularity



See also <https://www.universal-robots.com/articles/ur/application-installation/what-is-a-singularity/>

Resources



Find the latest UR software at [https://www.universal-robots.com/download/?filters\[\]=98763&query=](https://www.universal-robots.com/download/?filters[]=98763&query=) If experiencing problems with a robot, always try installing the latest software.

Manuals online: https://myur.universal-robots.com/manuals/content/SW_5_14/Documentation%20Menu

Release notes: <https://www.universal-robots.com/articles/ur/release-notes/release-note-software-version-514xx/>

Support articles: <https://www.universal-robots.com/articles/?>

Forum: <https://forum.universal-robots.com/>

Discord: <https://discord.com/invite/sEjRgEf6fp>

Inspect logs with Log Viewer: <https://www.universal-robots.com/download/software-e-series/support/ur-log-viewer/log-viewer-v1210/>

Collect realtimedata with
https://github.com/UniversalRobots/RTDE_Python_Client_Library

GitHub <https://github.com/orgs/UniversalRobots/repositories>

A young boy with blonde hair, wearing a white lab coat, is focused on a blue handheld device. He is standing in a large, brightly lit industrial space, likely a factory or exhibition hall. To his left, a large, grey and blue Universal Robots robotic arm is visible. In the background, several other people are standing and observing, suggesting a public demonstration or event. The floor is marked with yellow and black safety lines. The overall atmosphere is one of curiosity and engagement with technology.

Questions?

Thank you

Let's change the world!

